



Philips HS1/FRx Battery Expiration



Install Before date = Last day of October, 2021

- Philips M5070A battery has a 5-year **shelf life** from date of manufacture (date on battery).
- From date of installation (in the AED), this date becomes a typical 4-year **standby life** during which time the device performs daily, weekly and monthly self-tests (provided one installation, no use).
- The 4-year standby life applies for any battery installed prior to the “install before” date indicated on it.
- Below are three examples of when a battery like 2021-10 above could be installed*:



**NOTE: Before battery is totally depleted, the AED will start chirping when it reaches a low battery condition. Press the flashing i-Button to confirm the low battery condition, and replace with a fresh battery. Typically, this would be 4 years after battery initial installation. If you have not done so already, equip your AED with a spare battery today to avoid having to purchase a fresh one in a hurry.*